

# Sequencing Protocol for BIC Genomics Core Facility

## I. Reagents

BigDye Terminator v.3.1 or v.1.1 Cycle Sequencing Kit

## II. Storage

1. Store kit at -15 to -25°C.
2. Before each use, allow to thaw at room temp (do not heat).
3. Materials should be kept on ice during use.
4. Avoid repeated freeze/thaw; instead, aliquot into multiple tubes.

## III. Template Preparation Methods

1. Starting concentration of templates:
  - a. double-stranded plasmid less than 7K: 0.300µg-0.800µg
  - b. double-stranded plasmid 7K-12K: 0.800µg-1.500µg\*

\*Certain expression vectors sequence very poorly regardless of the concentration.

2. Genomic questions: see Zhao Li, 295-3657
3. Expected amount of PCR product:
  - a. 50bp-200bp: 10ng-30ng
  - b. 200bp-700bp: 30ng-60ng
  - c. 700bp-1000bp:60ng-100ng

## IV. Reaction Conditions

1. X µl of template
2. Y µl of H<sub>2</sub>O
3. 10 pmol primer (If concentration is unknown, more is better than less.)
4. 4.0µl of BigDye v.3.1 or v.1.1
5. Total volume=20µl
6. Mix well and spin briefly.

## V. Thermal Cycler Profile

Hot start; 25 cycles; 3-5min final extension at 60°C

Step	Action
1	Place tubes in thermal cycler and set the volume to 20µl.
2	Repeat the following for 25 cycles: <ol style="list-style-type: none"><li>a. Rapid thermal ramp to 96°C</li><li>b. 96°C for 10sec</li><li>c. Rapid thermal ramp to 50°C</li><li>d. 50°C for 5sec*</li><li>e. Rapid thermal ramp to 60°C</li><li>f. 60°C for 4min</li></ol>
3	Rapid thermal ramp to 4°C and hold until ready to purify.
4	Spin down the contents of the tubes in microcentrifuge.

\*Set temperature 1-2°C below T<sub>m</sub> of primer used.

## VI. Purification

Purify PCR product using a DTR gel filtration column.